




Neo SinusAll Kit



SCA Kit


SLA Kit

Surgical Kit

The SinusAll Kit combines NeoBiotech's SCA (Sinus Crestal Approach) Kit, SLA (Sinus Lateral Approach) Kit, and the Surgical Kit into a single kit. The product allows maxillary sinus implant procedures to take place using a single kit, rather than having to prepare a sinus kit and a surgical kit separately.

Crestal Approach Tool Components

• S-Reamer




Ø2.8 Ø3.2 Ø3.6 Ø3.9

SASR28 SASR32 SASR36 SASR39

Diameter

Product name

• Rotary elevator



Ø2.8 Ø3.2 Ø3.6 Ø3.9

SARE28 SARE32 SARE36 SARE39


• Stopper



Length 1mm 2mm 3mm 4mm 5mm 6mm 7mm 8mm 9mm 10mm 11mm 12mm


Product name CDS010 CDS020 CDS030 CDS040 CDS050 CDS060 CDS070 CDS080 CDS090 CDS100 CDS110 CDS120

- Depth gauge



Product name SADG01

- Bone carrier



Product name JUMB02

- Bone condenser




Product name SABCO1

- Aqua system (Aqua Tap)



Product name AQT

- Aqua system (Aqua valve, hose, syringe)



Product name AQST

- Drill Handle



Product name SADH01

Lateral Approach Tool Components

- Sinus Elevator



Product name SLE01



Product name SLE02



Product name SLE03

- Side Cutter



Product name STC30

- LS-Reamer



Product name

SALS640

- C-Guide Reamer & C-Reamer



Product name

SACG610



Product name

SAC640


- Lateral Stopper



Length	1mm	2mm	3mm	4mm
Product name	CDS010	CDS020	CDS030	CDS040

Implant Placement Tool Components

- Hex Driver



Product name HD1215S

- Fixture Driver




Length	Ratchet (Short)	Ratchet (Long)	Contra Angle
Product name	ISFD10R	ISFD15R	ISFD05C

- Torque wrench



Product name TW60

- Sinus Drill



Diameter	Ø2.2	Ø2.8	Ø3.2	Ø3.7	Ø4.0
Product name	SAD22	SAD28	SAD32	SAD37	SAD40